

Eppendorf Certificate

Multipette® E3/E3x and Repeater® E3/E3x

Multipette E3/E3x and Repeater E3/E3x are the enhanced successor models of the Multipette stream/Xstream and Repeater stream/Xstream. Multipette E3/E3x and Repeater E3/E3x do not differ in specification for any of the properties listed below:

- Technical specifications (accuracy and precision)
- Volume range from 1 µL to 50 mL
- Recognition of Combitips size and display of dispensing volume
- Possibility to dispense up to 100 times with one tip filling
- Positive displacement principle for dispensing
- Contamination free dispensing of hazardous liquids
- Precise and accurate dispensing of viscous solutions or liquids with high vapor pressure

The optimization of the Multipette E3/E3x and Repeater E3/E3x refers to the properties below.

- New sleep function for extending battery life
- Housing material with better chemical resistance
- New embedded RFID data chip with product properties
- Works only with improved Combitips advanced® for better tip attachment and more reliable wall dispensing
- Enlarged display size for a better readability
- Optimized tip sensor to detect a larger number of different tips and a better resistance against mechanical and chemical influences
- Save favorite settings
- New “Aspirate and Dispense” operation mode
- Individual setting of service intervals

Hamburg, July 19th, 2017



Brigitte Klose
Global Product Manager
Liquid Handling



Dr. Birgit Schreiber
Vice President
Quality Management & Regulatory Affairs

Your local distributor: www.eppendorf.com/contact
Eppendorf AG · 22331 Hamburg · Germany
E-mail: eppendorf@eppendorf.com

ISO 9001
Certified

ISO 13485
Certified

ISO 14001
Certified