

Safety Data Sheet

Trade name: EMKARATE

Current version : 1.0.0, issued: 12.01.2023

Replaced version: -, issued: -

Region: USA

SECTION 1: Identification

1.1 Product identifier

Trade name

EMKARATE

1.2 Recommended use of the chemical and restrictions on use

Recommended use of the chemical and restrictions on use

Lubricant for compressors for professional use in ultra-low temperature freezers.

Restrictions on use

No data available.

1.3 Details of the supplier of the safety data sheet

Address

Eppendorf North America, Inc.
175 Freshwater Blvd

Enfield, CT 06082
United States of America

Telephone no. +1 800 645 3050

Fax no. +1 516 334 7506

Information provided by / telephone

email: certificates@eppendorf,.com

Supplier

Address

Eppendorf SE
Barkhausenweg 1
D-22339 Hamburg
Germany

1.4 Emergency phone number

For medical and spill advice:

Phone: 800-451-8346 (Verisk 3E; client ID number: 12919)

SECTION 2: Hazard(s) identification

2.1 Classification of the chemical

Classification information

This product is not subject to classification in accordance with paragraph (d) of 29 CFR § 1910.1200 "Hazard Communication".

2.2 Label elements

Not relevant

2.3 Labelling information

No data available

2.4 Hazards not otherwise classified (HNOC)

No data available

SECTION 3: Composition / information on ingredients

3.1 Substances

Not applicable. The product is not a substance.

3.2 Mixtures

Chemical characterization

Mixture (preparation)

Safety Data Sheet

Trade name: EMKARATE

Current version : 1.0.0, issued: 12.01.2023

Replaced version: -, issued: -

Region: USA

Hazardous ingredients

This product does not contain substances to be mentioned in accordance with paragraph (d) of 29 CFR § 1910.1200 "Hazard Communication".

SECTION 4: First-aid measures

4.1 Description of necessary first-aid measures

General information

Remove contaminated clothing and shoes immediately, and launder thoroughly before reusing.

After inhalation

Remove affected person from the immediate area. Ensure supply of fresh air.

After skin contact

Wash off with soap and water. Get medical attention if pain still persists.

After eye contact

Remove contact lenses. Rinse eye thoroughly under running water keeping eyelids wide open and protecting the unaffected eye (at least 10 to 15 minutes). Get immediate ophthalmic treatment.

After ingestion

Do not induce vomiting - aspiration hazard. Never give anything by mouth to an unconscious person. If individual is drowsy or unconscious, place in recovery position (on left side, with head down). By continuous complaints consult a physician.

4.2 Most important symptoms and effects, both acute and delayed

No data available.

4.3 Indication of any immediate medical attention and special treatment needed

Treat symptomatically

SECTION 5: Fire-fighting measures

5.1 Extinguishing media

Suitable extinguishing media

Carbon dioxide; Dry chemical extinguisher; Water spray jet; Water mist

Unsuitable extinguishing media

Full water jet

5.2 Specific hazards arising from the chemical

Product floats on water.

5.3 Special protective equipment and precautions for fire-fighters

As in any fire, wear self-contained breathing apparatus (MSHA/ NIOSH approved or equivalent) and full protective gear.

SECTION 6: Accidental release measures

6.1 Personal precautions, protective equipment and emergency procedures

For non-emergency personnel

Refer to protective measures listed in sections 7 and 8. Ensure adequate ventilation.

For emergency responders

No data available. Personal protective equipment (PPE) - see Section 8.

6.2 Methods and materials for containment and cleaning up

Pick up with absorbent material (e.g., sand, sawdust, general-purpose binder). Send in suitable containers for recovery or disposal.

SECTION 7: Handling and storage

7.1 Precautions for safe handling

Precautions for safe handling

Ensure adequate ventilation.

Safety Data Sheet

Trade name: EMKARATE

Current version : 1.0.0, issued: 12.01.2023

Replaced version: -, issued: -

Region: USA

General protective and hygiene measures

Do not eat, drink or smoke during work time. Keep away from foodstuffs and beverages. Wash soiled clothing. Wash hands before breaks and after work.

Advice on protection against fire and explosion

No special measures necessary.

7.2 Conditions for safe storage, including any incompatibilities

Technical measures and storage conditions

Keep container tightly closed in a cool, well-ventilated place.

Requirements for storage rooms and vessels

Containers which are opened must be carefully resealed and kept upright to prevent leakage. Always keep in containers of same material as the original.

Incompatible products

Do not store together with: Oxidizing agents; Acids; Bases

SECTION 8: Exposure controls / personal protection

8.1 Control parameters

No parameters available for monitoring.

8.2 Exposure controls

Appropriate engineering controls

No data available.

Individual protection measures, such as personal protective equipment

Respiratory protection

In case of aerosol and mist formation, take appropriate measures for breathing protection in the event workplace threshold values are not specified.

Eye / face protection

Safety glasses recommended

Hand protection

Sufficient protection is given wearing suitable protective gloves in the event of risk of skin contact with the product. Before use, the protective glove should be tested in any case for its specific work-station suitability (i.e. mechanical resistance, product compatibility and antistatic properties). Adhere to the manufacturer's instructions and information relating to the use, storage, care and replacement of protective gloves. Protective gloves shall be replaced immediately when physically damaged or worn. Design operations thus to avoid permanent use of protective gloves.

Appropriate Material neoprene

Other protection measures

work clothing

SECTION 9: Physical and chemical properties

Form/Color

Liquid

light yellow

Odor

mild

Odor threshold

No data available

pH value

No data available

Boiling point / boiling range

Value > 250 °C

Melting point / melting range

No data available

Decomposition point / decomposition range

Safety Data Sheet

Trade name: EMKARATE

Current version : 1.0.0, issued: 12.01.2023

Replaced version: -, issued: -

Region: USA

No data available		
Pour point		
Value	-52	°C
Flash point		
Value	> 240	°C
Method	Cleveland open cup	
Auto-ignition temperature		
Value	390	°C
Oxidising properties		
No data available		
Explosive properties		
No data available		
Flammability		
No data available		
Lower flammability or explosive limits		
No data available		
Upper flammability or explosive limits		
No data available		
Vapor pressure		
Value	0.1	Torr
Reference temperature	20	°C
Vapor density		
No data available		
Evaporation rate		
No data available		
Relative density		
Value	0.995	
Reference temperature	20	°C
Density		
No data available		
Solubility in water		
Remarks	insoluble	
Solubility(ies)		
No data available		
Partition coefficient: n-octanol/water		
No data available		
Viscosity		
Value	18.9	mm²/s
Reference temperature	40	°C
Other information		
Other information		
No data available		

SECTION 10: Stability and reactivity

10.1 Reactivity

No data available.

10.2 Chemical stability

Stable under recommended storage and handling conditions(See section 7).

Safety Data Sheet

Trade name: EMKARATE

Current version : 1.0.0, issued: 12.01.2023

Replaced version: -, issued: -

Region: USA

10.3 Possibility of hazardous reactions

Dangerous reactions are not to be expected when handling product according to its intended use.

10.4 Conditions to avoid

None known.

10.5 Incompatible materials

Oxidizing agents; strong bases; strong acids

10.6 Hazardous decomposition products

Carbon monoxide and carbon dioxide

SECTION 11: Toxicological information

Information on toxicological effects

Acute oral toxicity

No data available

Acute dermal toxicity

No data available

Acute inhalational toxicity

No data available

Skin corrosion/irritation

No data available

Serious eye damage/irritation

No data available

Respiratory or skin sensitisation

No data available

Germ cell mutagenicity

No data available

Reproduction toxicity

No data available

Carcinogenicity

No data available

NTP Report on Carcinogens (RoC)

According to the data available and/or specifications supplied by upstream suppliers, this product does not contain any substance(s) listed "known" (to be a human carcinogen) or "RAHC" (reasonably anticipated to be a human carcinogen).

OSHA Standards listing cancer as a health effect (Table VII.6.2)

According to the data available and/or specifications supplied by upstream suppliers, this product does not contain any substance(s) listed as carcinogenic.

IARC Monographs on the Evaluation of Carcinogenic Risks to Humans

According to the data available and/or specifications supplied by upstream suppliers, this product does not contain any substance(s) classified by the IARC Monographs.

STOT - single exposure

No data available

STOT - repeated exposure

No data available

Aspiration hazard

No data available

Delayed and immediate effects as well as chronic effects from short and long-term exposure

Product specific toxicological data are not available.

Safety Data Sheet

Trade name: EMKARATE

Current version : 1.0.0, issued: 12.01.2023

Replaced version: -, issued: -

Region: USA

SECTION 12: Ecological information

12.1 Toxicity

Toxicity to fish (acute)
No data available
Toxicity to fish (chronic)
No data available
Toxicity to Daphnia (acute)
No data available
Toxicity to Daphnia (chronic)
No data available
Toxicity to algae (acute)
No data available
Toxicity to algae (chronic)
No data available

12.2 Persistence and degradability

Biodegradability			
No	Product name		
1	EMKARATE		
Value		76	%
Duration		28	day(s)
Method	OECD 301 B		
Evaluation	readily biodegradable		
Value			%

12.3 Bioaccumulative potential

No data available.

12.4 Mobility in soil

No data available.

12.5 Other adverse effects

No data available.

12.6 Other information

Other information
Do not discharge into drains or waters and do not dispose of in public landfills.

SECTION 13: Disposal considerations

13.1 Disposal considerations

Product

Dispose in accordance with federal, state and local regulations.

Packaging

Residuals must be removed from packaging and when emptied completely disposed of in accordance with the regulations for waste removal. Incompletely emptied packaging must be disposed of in accordance with federal, state and local regulations.

SECTION 14: Transport information

14.1 US DOT (Department of Transportation)

The product is not subject to DOT regulations.

14.2 Transport IMDG

The product is not subject to IMDG regulations.

14.3 Transport ICAO-TI / IATA

Safety Data Sheet

Trade name: EMKARATE

Current version : 1.0.0, issued: 12.01.2023

Replaced version: -, issued: -

Region: USA

The product is not subject to ICAO-TI / IATA regulations.

14.4 Other information

No data available.

14.5 Environmental hazards

Information on environmental hazards, if relevant, pls. see 14.1 - 14.3.

14.6 Special precautions for user

No data available.

14.7 Transport in bulk according to Annex II of MARPOL 73/78 and the IBC Code

Not relevant

SECTION 15: Regulatory information

15.1 US Federal Regulations

Toxic Substances Control Act (TSCA)

According to the data available and/or specifications supplied by upstream suppliers, all components of this product are on the TSCA Inventory or are exempted from notification requirements.

CERCLA (Comprehensive Environmental Response, Compensation and Liability Act)

According to the data available and/or specifications supplied by upstream suppliers, this product does not contain any substance(s) that are considered being subject to Table 302.4 "List of Hazardous Substances and Reportable Quantities".

Clean Air Act (CAA): Hazardous Air Pollutants (HAPS)

According to the data available and/or specifications supplied by upstream suppliers, this product does not contain any substance(s) that are considered being "Hazardous Air Pollutants".

Clean Water Act (CWA) Section 307 (Toxic Pollutants) and Section 311 (Hazardous Substances)

According to the data available and/or specifications supplied by upstream suppliers, this product does not contain any substance(s) that are considered to be "Hazardous Substances".

EPCRA SARA Section 313, Toxic Chemical Release Inventory

According to the data available and/or specifications supplied by upstream suppliers, this product does not contain any substance(s) that are considered being subject to "SARA Reporting".

Emergency Planning & Community Right-to-Know Act (EPCRA) - Sections 302 and 304

According to the data available and/or specifications supplied by upstream suppliers, this product does not contain any substance(s) that are considered being subject to reporting requirements.

Clean Air Act (CAA) - Section 602 Ozone depleting substances

According to the data available and/or specifications supplied by upstream suppliers, this product does not contain any substance(s) listed as ozone depleting substance.

15.2 US State Regulations

Pennsylvania State Right to Know Hazardous Substance List

According to the data available and/or specifications supplied by upstream suppliers, this product does not contain any substance(s) that are considered to be "Hazardous Substances".

New Jersey State Right to Know Hazardous Substance List

According to the data available and/or specifications supplied by upstream suppliers, this product does not contain any substance(s) that are considered to be "Hazardous Substances".

Massachusetts State Right to Know Hazardous Substance List

According to the data available and/or specifications supplied by upstream suppliers, this product does not contain any substance(s) that are considered to be "Hazardous Substances".

California Code of Regulations Title 8 §339, Hazardous Substances List

According to the data available and/or specifications supplied by upstream suppliers, this product does not contain any substance(s) that are considered to be "Hazardous Substances".

Right to Know - California State - Proposition 65

According to the data available and/or specifications supplied by upstream suppliers, this product does not contain any substance(s) that are considered to be "Hazardous Substances".

Safety Data Sheet

Trade name: EMKARATE

Current version : 1.0.0, issued: 12.01.2023

Replaced version: -, issued: -

Region: USA

SECTION 16: Other information, including date of preparation or last revision

16.1 Other information

No data available.

16.2 Date of preparation/revision of the SDS

Current version	1.0.0	12.01.2023
Replaced version	-	-

16.3 Department issuing safety data sheet

UMCO GmbH
Georg-Wilhelm-Str. 187, D-21107 Hamburg, Germany
Tel.: +49 40 / 555 546 300 Fax: +49 40 / 555 546 357 e-mail: umco@umco.de

16.4 Conclusions

This information is based on our present knowledge and experience.

The safety data sheet describes products with a view to safety requirements.

It does not, however, constitute a guarantee for any specific product properties and shall not establish a legally valid contractual relationship.

Alterations/supplements:

Alterations to the previous edition are marked in the left-hand margin.

Document protected by copyright. Alterations or reproductions require the express written permission of UMCO GmbH.

Prod-ID 779098

16.5 Key / Legend

Safety Data Sheet

Trade name: EMKARATE

Current version : 1.0.0, issued: 12.01.2023

Replaced version: -, issued: -

Region: USA

5E - indicates a substance that is the subject of a TSCA section 5(e) order
5F - indicates a substance that is the subject of a TSCA section 5(f) rule
ACGIH - American Conference of Governmental Industrial Hygienists
ATE – Acute Toxicity Estimate
BOD - Biochemical Oxygen Demand
CA/MA/MN/NJ/PA - California/Massachusetts/Minnesota/New Jersey/Pennsylvania*;
CAS - Chemical Abstracts Service
CERCLA - Comprehensive Environmental Response, Compensation, and Liability Act
CFR - Code of Federal Regulations (US)
CLP - Classification, Labelling, and Packaging
DOT - Department of Transportation Commission
ECHA – European Chemicals Agency
EC / EL – Effective Concentration / Effective Loading
EPA - Environmental Protection Agency
IARC - International Agency for Research on Cancer
IATA - International Air Transport Association
IBC – Intermediate Bulk Container
ICAO - International Civil Aviation Organization
IDL - Ingredient Disclosure List
IDLH - Immediately Dangerous to Life and Health
IMDG - International Maritime Dangerous Goods
IUCIID - International Uniform Chemical Information Database
Kow - Octanol/water partition coefficient
LD/LC - Lethal Dose/ Lethal Concentration
MARPOL – Marine Pollution
MEL - Maximum Exposure Limits
NFPA - National Fire Protection Agency
NIOSH - National Institute for Occupational Safety and Health
NJTSR - New Jersey Trade Secret Registry
NOEL/NOEC – No Observed Effect Level / No Observed Effect Concentration
NTP - National Toxicology Program
OSHA - Occupational Safety and Health Administration
PEL- Permissible Exposure Limit
PMN - indicates a commenced Pre Manufacture Notice substance
RCRA - Resource Conservation and Recovery Act
S - indicates a substance that is identified in a final Significant New Use Rule
SARA - Superfund Amendments and Reauthorization Act
SP - indicates a substance that is identified in a proposed Significant New Use Rule
STEL - Short-term Exposure Limit
STOT – Specific Target Organ Toxicity
T - indicates a substance that is the subject of a final TSCA section 4 test rule or order
TDG - Transportation of Dangerous Goods
TLV - Threshold Limit Value
TP - indicates a substance that is the subject of a proposed TSCA section 4 test rule or order
TSCA - Toxic Substances Control Act
TWA - Time Weighted Average
UN/NA - United Nations /North American
US - United States
VLE - Exposure Limit Value (Mexico)
VOC – Volatile Organic Compounds
XU - indicates a substance exempt from reporting under the Chemical Data Reporting Rule